

## + InSpec™ HA

Sterile Sporicidal Preparation of Hypochlorous Acid (HC10 $\delta$ OCI $\delta$  + H $\delta$ <sup>+</sup>). Low odour, safe handling, low residue and fully effective.

### + INFORMATION

#### Microbiological Minimum Efficacy Testing Contact

EN 1276 Bacteria	5 minutes	EN 13697 Bacteria	5 minutes
EN 1650 Fungi	15 minutes	EN 13697 Fungi	15 minutes
EN 13704	1 minute (Bacillus Subtilis)		



### + ORDERING INFORMATION

#### InSpec™ HA Products

Product	Pack Size	Case Size	Redditch Code	AMD Code
InSpec™ HA 1L Trigger	1L Trigger Spray	6 x 1L	HAPW50-1LS	D06C10H
InSpec™ HA 5L Screw Cap	5 Litre Screw Cap	2 x 5L	HAPW50-5LS	D02C00H
InSpec™ HA Burstable Pouch	15 wipes/pouch	20 pouches	HAPW50BP15S	W10C20H

Please use the AMD Code when ordering.

AMD are the official UK distributor for Redditch Medical's range of products to the NHS, pharmacies and compounders.

### + PRODUCT INFORMATION

#### Use Instructions

InSpec HA™ is designed for spraying, wiping and mopping applications. Hold spray approximately 15cm to 20cm from area to be treated. 5L versions, pour into an appropriate container for mopping.

#### Technical Data

Appearance: Colourless aqueous solution with a characteristic mild chlorine odour. pH 4.0 to 6.5.

#### Storage

Store upright in original closed containers away from sunlight and extremes of temperature. Full guidance on the handling and disposal of this product is available in the Safety Data Sheet (SDS).

#### Safe Handling

Always wear gloves and goggles or face protection. Always read the label and SDS before use.

#### Material Compatibility

Application of solutions, when used as directed, will not affect materials normally encountered in the cleanroom. See compatibility information in the Technical File.

#### Benefits

- One of the safest and fully effective sporicides. All formats are non hazardous
- Very low residue, corrosion and odour
- Stabilised solution - 12 months shelf life
- 6 months validated 'in-use' shelf life
- Contents protected using 'bag-in-bottle' presentation with premium range (spray bottle)
- 1L and 5L formats available as Ready to Use
- Multiple bags for cleanroom transfer
- Ideal rotational partner for InSpec™ QT, InSpec™ QT+ and InSpec™ AN
- Manufactured to GMP
- Active Substance compliant with Article 95 of the Biocidal Products Regulation (BPR)

#### Technical User Information

InSpec™ HA is a sterile sporicide available ready to use.

InSpec™ HA is manufactured to GMP and is filled and capped in an ISO class 5 environment.

The solution is filtered through a 0.2 micron filter and aseptically filled into pre-gamma irradiated packaging (25-45kGy) and tested to give a sterility assurance level SAL of 10<sup>-6</sup>.

For screw cap bottles use the entire contents after opening.

#### Aseptic Medical Devices

Newmarket Drive, Derby DE24 8SW | +44 (0) 1332 755622  
sales@asepticmedical.com | www.asepticmedical.com

Aseptic Medical Devices is a trading name of Riverside Medical Packaging Company Limited