

AMD Sterile Transfer Sets & Accessories


# Terumo Needles for Syringe Filling Formatted for Aseptic Pharmacy




## Overview




### How to access the hub of the K-Pack Terumo Needle



1. Grip needle case at both ends



2. Make a slight breaking movement (at the perforation) on both ends



3. Pull the cap from the needle

AMD selected the Terumo K-Pack II needle for aseptic transfer of liquids, as it is specifically designed to allow insertion into a secondary packaging.

Each needle is individually encapsulated in an ergonomically shaped, colour-coded hard cap and case. This format ensures safety and compliance with good aseptic technique.

### Ordering Information

  Stock Item      Made to Order

Needles - Pouched & Double Outer Wrapped		Quantity per Pouch				
		5	10	25	50	100
18 Gauge	Encapsulated Needle - 38mm (1.5") in Length	N05A18T	N10A18T	N25A18T	N50A18T	N00A18T
19 Gauge	Encapsulated Needle - 38mm (1.5") in Length	N05A19T	N10A19T	N25A19T	N50A19T	N00A19T
21 Gauge	Encapsulated Needle - 38mm (1.5") in Length	N05A21T	N10A21T	N25A21T	N50A21T	POA
23 Gauge	Encapsulated Needle - 32mm (1.25") in Length	N05A23T	N10A23T	N25A23T	N50A23T	POA