



+ AMD Isolator Cleaning Tool and Bonnet

Formatted for Aseptic Pharmacy

AMD Sterile Transfer Sets & Accessories



+ ORDERING INFORMATION

AMD Isolator Cleaning Tool and Bonnet	
Description	Product Code
AMD Isolator Cleaning Bonnet, Individually Pouched (carton of 25)	X70A25R
AMD Isolator Cleaning Bonnet, Individually Pouched (carton of 50)	X70A50R
AMD Isolator Cleaning Tool	X70B01R

+ PRODUCT INFORMATION

AMD's Isolator Cleaning Tool and Cleaning Bonnets are specifically designed for isolator, laminar flow cabinet cleaning and other cleanroom areas where access is restricted.

+ CLEANING TOOL

The triangular head shape of the cleaning tool allows for improved access in tight corners ensuring areas are cleaned effectively and efficiently.

- Fully autoclavable
- Made from 316 stainless steel with aluminium handle for optimum reach
- Easy to use, lightweight and manoeuvrable

+ CLEANING BONNET

Designed to be used with the cleaning tool for efficient cleaning of environments that are critical to the quality of compounded products. The presentation of a dry absorbent bonnet allows for customers to use their preferred biological agents and cleaning solutions.

- Single use disposable low particulate bonnets are manufactured from a 2 ply polyester blend material
- Elastic sewn edge enables a secure fit on the isolator cleaning tool head during use
- Supplied triple sterile wrapped and packed within a class 7 cleanroom sterilised by gamma irradiation
- AMD's unique packaging ensures easy opening, even when wearing gloves

Aseptic Medical Devices

Newmarket Drive, Derby DE24 8SW | +44 (0) 1332 755622
 sales@asepticmedical.com | www.asepticmedical.com

Aseptic Medical Devices is a trading name of Riverside Medical Packaging Company Limited

