



+ Individual Sachet IPA 70% Wipe

Formatted for Aseptic Pharmacy
AMD Sterile Transfer Sets & Accessories



+ ORDERING INFORMATION


Individual Sachet IPA 70% Wipe		Quantity per Pouch		
Description	Size	5	10	25
AMD Sterile Wipe with 70% IPA	200mm x 120mm	W05A02A	W10A02A	W25A02A

+ PRODUCT INFORMATION

In the course of aseptic manufacturing, there are small surfaces which require cleaning. To enable this activity in LAF cabinets and isolators AMD's range of small sterile wipes has been developed to meet this need.

+ B Braun Wipe

Formatted for Aseptic Pharmacy
AMD Sterile Transfer Sets & Accessories



+ ORDERING INFORMATION

Individual Sachet IPA 70% Wipe		Quantity per Pouch	
Description	Size	5	10
B Braun Alcohol Wipes with 70% IPA	30mm x 65mm	W05A03B	W10A03B

+ PRODUCT INFORMATION

B Braun Alcohol Pads are 30mm x 65mm and folded in half to fit the individual sealed foil pouch. The non-linting wipe is impregnated with 70% IPA.

Aseptic Medical Devices

Newmarket Drive, Derby DE24 8SW | +44 (0) 1332 755622
sales@asepticmedical.com | www.asepticmedical.com

Aseptic Medical Devices is a trading name of Riverside Medical Packaging Company Limited

