




## + Waste Bags

Formatted for Aseptic Pharmacy  
**AMD Sterile Transfer Sets & Accessories**



## + ORDERING INFORMATION

Waste Bags	Quantity per Pouch	
	10	20
<b>Description</b>		
Sterile Polyethylene Zip Lock Bag 152mm x 228mm	Z10B05W	POA
Sterile Polyethylene Zip Lock Bag 210mm x 297mm	Z10B10W	POA
Sterile Polyethylene Zip Lock Bag 300mm x 405mm	Z10B20W	Z20B20W
Large Transfer Waste Bags 600mm x 900mm	X10A30A	POA

## + PRODUCT INFORMATION

In the course of aseptic manufacture a significant non-sharps waste stream can be generated. For ease of removal from LAF cabinets or VHP isolators, AMD has developed a range of sterile polyethylene 'Zip Lock' waste bags and also provide a much larger sterile (non zip lock) bag for use in outer preparation areas.

These products are carefully folded and placed into triple packing using Aseptic Medical Devices' own specialist pouch: PeroxyMet/PeroxyClear.

### Aseptic Medical Devices

Newmarket Drive, Derby DE24 8SW | +44 (0) 1332 755622  
 sales@asepticmedical.com | www.asepticmedical.com

Aseptic Medical Devices is a trading name of Riverside Medical Packaging Company Limited

