

+ InSpec™ DE

Sterile Wipes with Denatured Ethanol 70% v/v and WFI Quality Water. Low residue and low linting, rapid action, fast drying, disinfectant.

+ INFORMATION

Microbiological Minimum Efficacy Testing Contact

EN 1276 Bacteria	5 minutes	EN 13697 Bacteria	5 minutes
EN 1650 Fungi	15 minutes	EN 13697 Fungi	15 minutes

+ ORDERING INFORMATION

InSpec™ DE Products

Product	Pack Size	Case Size	Redditch Code	AMD Code
InSpec™ DE Pouch 15pc	15 Wipes	25 per case	DEWFIPCH15S	W25C15T
InSpec™ DE Pouch 50pc	50 Wipes	16 per case	DEWFIPCH50S	W15C51T
InSpec™ DE Standard Tub	200 Wipes	12 per case	DEWFI12TB200S	W12C21T

Please use the AMD Code when ordering.

AMD are the official UK distributor for Redditch Medical's range to pharmacies and compounders.



+ PRODUCT INFORMATION

Use Instructions

InSpec™ DE Wipes are designed for wiping applications. Apply to surface to ensure complete coverage.

Storage

Store upright in original closed containers away from sunlight and extremes of temperature. Full guidance on the handling and disposal of this product is available in the Safety Data Sheet (SDS).

Safe Handling

Always wear gloves and goggles or face protection. Always read the label and SDS before use.

Material Compatibility

Application of solutions, when used as directed, will not affect materials normally encountered in the cleanroom. See compatibility information in the Technical File.

Pouch Wipes

InSpec™ DE Pouch Wipes are a high specification, low linting cleanroom wipes in a pouch.

DE Standard Pouch Wipes

Available as 15 or 50 wipes per pouch (55% cellulose; 45% polyester). Wipe dimension: 230mm x 230mm (9" x 9").

Tub Wipes

InSpec™ DE Tub Wipes, an economical solution for wiping applications.

DE Standard Tub Wipes

Wipes per tub: 200 x 23g cellulose. Dimensions of the wipe: 200mm x 200mm.

DE Premium Tub Wipes

Wipes per tub: 130 x 45g Polycellulose. Dimensions of the wipe: 200mm x 185mm.

Technical User Information

InSpec™ DE wipes are sterile ready to use disinfectant. This is a 70% v/v alcohol solution in WFI quality water.

InSpec™ DE wipes are manufactured to GMP in a cleanroom.

The solution is filtered through a 0.2 micron filter and gamma irradiated to a validated process (at 25-45kGy) to give a sterility assurance level SAL of 10⁻⁶.

InSpec™ DE is supplied with Certificates of Conformity, Analysis, Irradiation Sterility and Endotoxin.

All formats of InSpec™ DE have a stabilised 24 month shelf life.

Redditch Medical

Discovery 2, NETPark, 2 William Armstrong Way, Sedgefield, County Durham, TS21 3FD
UK Tel: +44 (0) 1332 755622 | International Tel: +44 (0) 1527 830942
sales@redditchmedical.com | www.redditchmedical.com