

## **+ Polyester Microfibre Mop with PU Foam**

High quality, efficient and effective



### **+ ORDERING INFORMATION**

#### **InSpec™ Impregnated Mops**

Product	Pack Unit	Redditch Code	AMD Code
InSpec™ IPA Impregnated Mop	2 Mops x 15	IPAWFIMOP-2S	MOPR02P
InSpec™ DE Impregnated Mop	2 Mops x 15	DEWFIMOP-2S	MOPR02T

#### **InSpec™ Burstable Mops**

Product	Pack Unit	Redditch Code	AMD Code
InSpec™ HA Burstable Mop	2 Mops x 10	HAPWMOP2S	MOPR02H
InSpec™ OX Burstable Mop	2 Mops x 10	OXPWMOP2S	MOPR02X
InSpec™ AN Burstable Mop	2 Mops x 10	ANWFIMOP2S	MOPR02A
InSpec™ QT Burstable Mop	2 Mops x 10	QTFWIMOP2S	MOPR02Q
InSpec™ QT+ Burstable Mop	2 Mops x 10	QTPLWFIMOP2S	MOPR02QT
InSpec™ N10 Burstable Mop	2 Mops x 10	N10WFIMOP2S	MOPR02N
InSpec™ N7 Burstable Mop	2 Mops x 10	N7WFIMOP2S	MOPR02D



Please use the AMD Code when ordering.

AMD are the official UK distributor for Redditch Medical's range to pharmacies and compounders.

### **+ PRODUCT INFORMATION**

#### **Description**

This Polyester-Microfibre Mop with felt medium to offer superb absorbency is designed for applying disinfectants to both horizontal and vertical surfaces, and is ideal for general cleaning and disinfection of ceiling, walls, and floors in the most critical environment. The microfibre material smoothly glides over floors and textured surfaces and provides superior contamination removal.

#### **Features**

- Knitted microfibre surface for good particle pick up
- Constructed with felt medium for high absorbency
- Ideal for cleaning and disinfection of Cleanrooms
- Chemical and abrasion resistant
- Light weight
- Low levels of particles
- Good gliding properties
- Lot to Lot traceable
- Sterile

#### **Applications**

- Excellent for spill control and cleaning large surfaces
- Compatible with most common detergents and disinfectants
- Suitable for use in Grade A-D cleanrooms

#### **Technical Data**

**Length:** 398 mm (±5 mm)

**Width:** 122 mm (±5 mm)

**Front material:** Polyester-Microfibre

**Backing material:** Nylon

**Absorbency capacity:** >305ml/mop (>6280ml/m²)

**Particles:** IEST-RPcc003.2 Helke Drumtest <3500 counts per mop >0.5µm

**Weight:** 900g (+/- 0.045)

#### **Redditch Medical**

Discovery 2, NETPark, 2 William Armstrong Way, Sedgefield, County Durham, TS21 3FD  
UK Tel: +44 (0) 1332 755622 | International Tel: +44 (0) 1527 830942  
sales@redditchmedical.com | www.redditchmedical.com