



+ Amber Syringe
Formatted for Aseptic Pharmacy
Sterile Transfer Sets & Accessories

Amber Syringe



Image shown for reference purposes only - product may differ in colour.

+ ORDERING INFORMATION

Amber Syringe		Quantity per Pouch			Certification	
Size	Description	5	10	25	CE	UKCA
50ml	Amber Luer Lock Syringe Fitted with Natural Coloured Stopper	S05A50X	S10A50X	S25A50X	✓	✓

+ PRODUCT INFORMATION

The Amber Syringe has a coloured barrel to reduce U.V. light transmission.

Testing has confirmed that the transmission does not exceed 15% at any wavelength between 290nm to 450nm.

The product is UKCA/CE marked and manufactured under ISO 13485 QMS.

AMD are always able to offer bespoke pack solutions such as different quantities per pouch or combinations of devices to meet specific procedure requirements bringing flexibility and scalability to any manufacturing process.

All syringe tips have a Luer Lock tip protector fitted to maintain the fluid pathway integrity.

The sterile transfer sets utilise Aseptic Medical Devices' own specialist pouch: PeroxyMet/PeroxyClear - a purpose designed pouch for aseptic area use - the pouch is alcohol, VHP and oxidising sporicide impermeable, yet transparent with a contrasting background for easy pre-use QC.

Riverside are the legal manufacturer of this product and take full responsibility for its safety and performance in accordance with the stated intended use.

Aseptic Medical Devices

Newmarket Drive, Derby DE24 8SW | +44 (0) 1332 755622
sales@asepticmedical.com | www.asepticmedical.com

Aseptic Medical Devices is a trading name of Riverside Medical Packaging Company Limited

